

P O W E R C H A I R



Seat Colors
Charcoal/Gray Vinyl



Candy Apple Red



Onyx



Body Colors
Blue
Hunter Green



**HD Version
450 lb.
weight capacity!**

Roll

golden
COMPASS™
HD

Stable platform, heavy duty center-wheel drive power chair with an articulating chassis for the ultimate in maneuverability and stability.



Rise. Roll. Rest. Ride.

GOLDEN
technologies

YOUR LIFE IN MOTION

Golden Compass™ HD Specifications

Model	GP620CC	GP620SS
Medicare Code	K0825	K0824
Weight Capacity	450 lbs	450 lbs
Drive Wheels	Center	Center
Maximum Speed	4 mph	4 mph
Operating Range ¹	12 Miles	12 Miles
Ground Clearance to transaxle/motor	3.5"	3.5"
Turning Radius	20"	20"
Type Batteries	2-NF22	2-NF22
Colors	Red, Blue, Green, Onyx	Red, Blue, Green, Onyx
Freewheel Mode	Yes	Yes
Electronic Speed Control	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes
Charger	Off Board	Off Board
Controller	Dynamic 75 Amp Remote	Dynamic 75 Amp Remote
Dimensions:		
Length	39.8"	39.8"
Width	25"	25"
Height (ground to top of back of seat)	38.25" - 42.75"	38.25" - 42.75"
Height (ground to top of back of headrest)	47.25" - 53.75"	47.25" - 53.75"
Ground to Top of Seat	21.75" - 24.25"	17.75" - 22.25"
Number of Seat Height Positions	5	5
Size of Increments	0.5"	0.5"
Weight of Unit: (assembled)		
Base	101 lbs	101 lbs
Seat w/Arms	51 lbs	48 lbs
Batteries	37 lbs ea	37 lbs ea
Tires:		
Rear	10" X 3.5"	10" X 3.5"
Casters	6" X 2"	6" X 2"
Standard Seat:		
Back Height - no headrest	Captain's High-back	Captain's Van/Pan
Back Height with headrest	21"	23"
Back Height with headrest	28.5" - 32.5"	30.5" - 32.5"
Width x Depth	20" X 20"	20" X 20", 22" X 20", 24" X 20"
Color	Coffee/Sand, Charcoal/Gray	Coffee/Sand, Charcoal/Gray
Footrest	Yes	Yes
Height adjustable	3" - 7"	3" - 7"
Angle adjustable	0 - 35 degrees	0 - 35 degrees
Warranty:		
Frame	5 Years	5 Years
Drive Train	1 Year	1 Year
Electronics ²	1 Year	1 Year

¹. Battery range will vary due to rider weight, drive surface and terrain.
². Electronics warranty excludes batteries.

4 exciting colors



Candy Apple Red



Blue



Hunter Green



Onyx



Power Chair Accessories and Options

- Captain's seat (GP620CC)
- Pro Flex seat (GP620SS)
- Van/pan seat (GP620SS)
- Swing-away legrests
- Elevating legrests
- Stump support
- Footplate extender
- Cane/crutch holder kit
- Walker holder kit
- Oxygen tank holder kit
- Pack N Go



Compass GP620SS shown