

Spring 2023

nova

Bathroom Safety



Explore all of NOVA's products at:
novajoy.com

Do You Know the BATHROOM HOT SPOTS?



We use the bathroom every day ... 365 days of each year, yet this room is often the most dangerous place in your home. Being proactive with effective product solutions is the key to fall prevention.

For less than a dollar a day, you can equip your bathroom with products to keep you safe. Make your bathroom happy and safe every day.

GET SOME LEVERAGE



HOT SPOT #1

If you have trouble getting on and off a chair, then you are also having trouble getting on and off a toilet. An elevated toilet seat and hand rails can make all the difference.

GETTING IN & OUT OF THE TUB



HOT SPOT #2

Getting in and out of the tub or shower can be tricky, and holding onto the sink or towel rack isn't safe. Using grab bars is a much better solution.

BATHING INDEPENDENTLY



HOT SPOT #3

Our tub or shower is our sanctuary, but without the proper safety equipment, this is another bathroom hot spot. Bath seats, hand-held shower heads and bath mats are a must!

VISIT US ON YouTube

NOVA's **YouTube** channel is jam packed with **great information, inspiration** and that **NOVAjoy**. We have **over 140 videos** and **counting** to educate, inform and inspire our customers.

Please do **Subscribe to our Channel** and feel free to give us a "thumbs up" or feedback. And, thank you for watching and sharing the joy!



Scan here to see our YouTube Channel!



GRAB BARS

Wall Grab Bars 8012-R to 8032-R



Chrome

Ivory

ITEM #	DESCRIPTION
8012CH-R	12" Wall Grab Bar - Chrome
8016CH-R	16" Wall Grab Bar - Chrome
8018CH-R	18" Wall Grab Bar - Chrome
8024CH-R	24" Wall Grab Bar - Chrome
8032CH-R	32" Wall Grab Bar - Chrome
8012PC-R	12" Wall Grab Bar - Ivory
8016PC-R	16" Wall Grab Bar - Ivory
8018PC-R	18" Wall Grab Bar - Ivory
8024PC-R	24" Wall Grab Bar - Ivory
8032PC-R	32" Wall Grab Bar - Ivory

- 1.5" wall stand-off
- 2.75" fixed flange
- 1" diameter tubing
- Knurled chrome grip (Chrome only)

Home Grab Bar 8609W-R

- Sleek, modern design with low profile
- Non-slip grip on bar promotes stability
- 250 lb weight capacity
- Only requires one wood stud to install and can be used throughout the home
- Best used in home danger zones such as bathrooms, doorways & stairwells
- Dimensions: 12.5"w x 2"d x 2.5"h
- Weight Capacity: 250 lb



Suction Cup Balance Assist Bar


8212-R

- Installs quickly; perfect for home or travel use
- Provides added support when entering and leaving a bath tub or shower
- Easy to attach and remove
- Secure suction indicator displays green for secure and red for not secure
- No tools required
- Dimensions: 5"w x 3.5"d x 12"h



16" Suction Cup Balance Assist Bar

8216B-R / 8216N-R



- Green and red indicators help ensure proper installation
- Waved tubing promotes secure hand grip
- Strong, easy-to-use suction cup mechanism
- Use for balance assist only, not for body weight leverage
- Each suction cup requires only 5" x 5" of smooth, non-porous tile space
- Dimensions: 4.25"w x 3.25"d x 20.5"h

Hand Held Shower Set

9301-R (2 Settings)

- Extra long 84" reinforced hose allows user to shower while seated
- On/Off switch allows for water flow control
- Includes wall mounting holder



9301-R

Hand Held Shower Sets

9311-R & 9312-R (5 Settings)

- Extra long 84" reinforced hose allows user to shower while seated
- On/Off switch allows for water flow control
- Includes wall mounting holder
- 5 Settings: Regular / Shower Spray / Massage / Water Drip / Shower Stop



9311-R
(White)

9312-R
(Chrome)

Tub Grab Bar

8220-R



- Provides added support when entering and leaving a bath tub or shower
- Made of 1" powder coated steel tubing
- Rubber contact points
- Clamp adjusts from 3" - 7"w
- Dimensions: 6"w x 11"d x 19"h
- Weight Capacity: 250 lb

Holder for Hand Held Shower

9306-R



- Attaches to most bath seat arms (3/4" diameter tube) or legs (1" diameter tube)
- Securely holds hand held shower in place for safety and convenience

Suction Cup Showerhead Holder

8213-R



- Showerhead is in easy reach when seated or bathing children
- Suction cup can be quickly removed and relocated as often as you like
- Convenient when bathing pets or cleaning out the tub
- Works on surfaces that are flat, smooth and non-porous such as tile, glass, mirror, fiberglass and acrylic
- Dimensions: 2.5"w x 4.25"d x 3.75"h

Shower Diverter Valve

9305-R



- Designed to divert water between your shower head and hand-held shower

Rubber Bath Mats


9350-R, 9351-R & 9351W-R

9350-R



- Attaches securely to tub surface and helps prevent slipping
- Anti-fungal treated
- 15.75"w x 27.5"l

9351-R



Also available in white (9351W-R)

Suction Cup Towel Hook

8214N-R (Nickel) & 8214B-R (Bronze)



8214N-R



8214B-R

- Easily add shower accessories to your bathroom
- Attaches via powerful suction mechanism without any tools
- Stylish finish for a modern touch to any bath or shower
- Easy to move or reposition

Bath Seat Cushion

9601W-R



- Adds comfort while bathing
- Non-slip Velcro strap prevents cushion from sliding
- Fits all standard sized bath seats
- Color: White
- Dimensions: 14.75"w x 8.5"d x 1"h

Bath Seat Basket

9605-R



- Holds bath essentials conveniently within reach to avoid bending, turning & twisting
- Easily attaches to the seat frame with included hardware
- Fits all standard size bath benches with 7/8" or 1" tube diameter
- Coated steel wire is rust resistant
- Dimensions: 12"w x 4.5"d x 4.25"h

Bath Seat
9039 with Arms & U-Shaped Cutout



- U-shaped cutout for improved hygiene
- Great for travel and easy to store
- Tool-free assembly in minutes
- Arms provide support for sitting & rising
- Adjustable height aluminum frame
- Latex-free rubber tips

• Seat Height Adjustment: 16" - 21"	• Leg Footprint: High - 19.5"w x 18.5"d Low - 18.75"w x 17.25"d
• Seat Dimensions: 15.75"w x 13"d	• Weight Capacity: 300 lb

Heavy Duty Bath Seat
9051 with U-Shaped Cutout



- U-shaped cutout for improved hygiene
- Adjustable height aluminum frame
- Latex-free rubber tips

• Seat Height Adjustment: 15.75" - 21.75"	• Leg Footprint: High - 20"w x 18.5"d Low - 18.75"w x 16.75"d
• Seat Dimensions: 19"w x 12.5"d	• Weight Capacity: 500 lb

Heavy Duty Bath Seat
9061 with Back & U-Shaped Cutout



- U-shaped cutout for improved hygiene
- Adjustable height aluminum frame
- Latex-free rubber tips
- Removable back

• Seat Height Adjustment: 15.75" - 21.75"	• Leg Footprint: High - 20"w x 18.5"d Low - 18.75"w x 16.75"d
• Seat Dimensions: 19"w x 12.5"d	• Weight Capacity: 500 lb


Bath & Shower Seats
9010 & 9020



- Built-in hand held shower holder
- Adjustable height aluminum frame
- Latex-free rubber tips
- Removable back (9020 only)

• Seat Height Adjustment: 16" - 22"	• Leg Footprint: High - 20"w x 17.75"d Low - 18.5"w x 16.5"d
• Seat Dimensions: 20"w x 12.25"d	• Weight Capacity: 300 lb

Bath & Shower Seats
9050 & 9060 with U-Shaped Cutout



- U-shaped design for cleaning purposes
- Adjustable height aluminum frame
- Latex-free rubber tips
- Removable back (9060 only)

• Seat Height Adjustment: 14.25" - 21.25"	• Leg Footprint: High - 20"w x 18"d Low - 18.5"w x 16.75"d
• Seat Dimensions: 19"w x 12.75"d	• Weight Capacity: 300 lb

Bath & Shower Seats
9016 & 9026



- Arms provide support for sitting and rising
- Adjustable height aluminum frame
- Latex-free rubber tips
- Removable back (9026 only)

• Seat Height: 16" - 20.75"	• Leg Footprint: High - 19.25"w x 18.25"d Low - 18.5"w x 17"d
• Seat Dimensions: 15.75"w x 13"d	• Weight Capacity: 300 lb
• Width Between Arms: 18.75"	

Deluxe Bath & Shower Seats
9120-R & 9100-R

- Safety handle for extra support and balance can be attached to either side.
- Built-in hand held shower holder
- Adjustable height aluminum frame
- Latex-free rubber tips



Red Safety Handle: Item # 9201-R

- | | |
|-----------------------------------------------|----------------------------------------------------------------------------|
| • Seat Height Adjustment:
15" - 21" | • Leg Footprint:
High - 17"w x 18.5"d
Low - 16.75"w x 17.25"d |
| • Seat Dimensions:
19.5"w x 16"d | • Weight Capacity: 300 lb |

Heavy Duty Bath & Shower Seats
9013 & 9023

- Adjustable height aluminum frame
- Built-in hand held shower holder
- Latex-free rubber tips
- Removable back (9023 only)



- | | |
|-----------------------------------------------|----------------------------------------------------------------------------|
| • Seat Height Adjustment:
17" - 23" | • Leg Footprint:
High - 19.75"w x 18"d
Low - 18.5"w x 16.75"d |
| • Seat Dimensions:
20"w x 12"d | • Weight Capacity: 500 lb |

Bath Seat
9131-R with Back

- Tool-free assembly in under 5 minutes
- Ergonomic, oversized seat supports up to 500 lb
- Stylish and modern colors give a classic feel to your bathroom
- Slip-resistant rubber tips add stability
- Non-corrosive aluminum legs will not rust



- | | |
|---------------------------------------------------|---------------------------------------------------------------------------|
| • Seat Height Adjustment:
16.5" - 20.5" | • Leg Footprint:
High - 19.25"w x 16.5"d
Low - 18.5"w x 16"d |
| • Seat Dimensions:
17"w x 14.75"d | • Weight Capacity: 500 lb |

Bath & Shower Seats
9121-R & 9101-R

- Great for travel and easy to store
- Tool-free assembly in minutes
- Built-in hand held shower holder
- Adjustable height aluminum frame
- Latex-free rubber tips



- | | |
|-----------------------------------------------|-------------------------------------------------------------------------|
| • Seat Height Adjustment:
17" - 22" | • Leg Footprint:
High - 20"w x 17.5"d
Low - 18"w x 16.75"d |
| • Seat Dimensions:
20"w x 12"d | • Weight Capacity: 300 lb |

Bath & Shower Seats
9035-R & 9036-R

- Great for travel and easy to store
- Tool-free assembly in minutes
- Arms provide support for sitting & rising
- Adjustable height aluminum frame
- Latex-free rubber tips



- | | |
|---------------------------------------------|-------------------------------------------------|
| • Seat Height: 16" - 21" | • Leg Footprint:
High - 19.75"w x 18" |
| • Seat Dimensions:
15.75"w x 13"d | • Low - 19"w x 16.75" |
| • Width Between Arms: 19" | • Weight Capacity: 300 lb |

Bath Seat
9037-R with Arms & U-Shaped Cutout

- Great for travel and easy to store
- U-shaped cutout for improved hygiene
- Tool-free assembly in minutes
- Arms provide support for sitting & standing
- Adjustable height aluminum frame
- Latex-free rubber tips



- | | |
|---------------------------------------------|------------------------------------------------------------------------------|
| • Seat Height: 16" - 21" | • Leg Footprint:
High - 19.5"w x 18.5"d
Low - 18.75"w x 17.25"d |
| • Seat Dimensions:
15.75"w x 13"d | • Weight Capacity: 300 lb |
| • Width Between Arms: 18.75" | |

Bath Stool
9006-R



- Great for travel and easy to store
- Tool-free assembly in minutes
- Adjustable height aluminum frame
- Latex-free rubber tips

• Seat Height Adjustment: 15.5" - 20.5"	• Leg Footprint: High - 15"w x 14.5"d Low - 14"w x 13.75"d
• Seat Dimensions: 12.5"w x 12.75"d	• Weight Capacity: 300 lb

Bath Stool
9005



- Great for use in small areas
- Adjustable height aluminum frame
- Latex-free rubber tips

• Seat Height Adjustment: 13.75" - 20"	• Leg Footprint: High - 15.5"w x 14.75"d Low - 14"w x 13.75"d
• Seat Dimensions: 12.75"w x 12.75"d	• Weight Capacity: 250 lb

Foldable Shower Chair
9400 with U-Shaped Cutout



- Great for travel/easy to store
- Folds and unfolds quickly
- U-shaped design for cleaning purposes
- Arms provide support for sitting and rising
- Adjustable height aluminum frame
- Latex-free rubber tips

• Seat Height Adjustment: 18" - 20.75"	• Leg Footprint: High - 21"w x 19.25"d Low - 20.5"w x 19"d
• Seat Dimensions: 16"w x 16"d	• Weight Capacity: 250 lb

Tub Shower Board
9200



- Width adjusts to attach securely to most bathtubs
- Built-in soap holder
- Rubber covered clamps resist slipping and protect surfaces from scratching
- Safety handle can be attached to either side for support

• Seat Dimensions: 31"w x 14.5"d	• Weight Capacity: 250 lb
• Width Adjustment: 17.75" - 28.5"	

Wall Mounted Shower Seat
9404



- Mounts directly to shower wall
- Saves space by folding up when not in use
- Includes wall mount anchors
- Adjustable height aluminum frame
- Latex-free rubber tips

• Seat Height Adjustment: 19.75" - 22.25"	• Overall Footprint: High - 15.75"w x 13.75"d Low - 15.75"w x 13.75"d
• Seat Dimensions: 16"w x 16"d	• Weight Capacity: 250 lb

Deluxe Transfer Bench
9070

- Built-in soap holder and hand held shower holder
- Converts easily for left or right hand transfers
- Adjustable height aluminum frame
- Latex-free rubber tips & suction tips



• Seat Height Adjustment: 18.75" - 23.25"	• Leg Footprint: High - 33"w x 17.5"d
• Seat Dimensions: 31.5"w x 14.5"d	• Weight Capacity: 300 lb

Standard Transfer Bench
9071

- Converts easily for left or right hand transfers
- Adjustable height aluminum frame
- Latex-free rubber tips & suction tips



• Seat Height Adjustment: 17.25" - 21.25"	• Leg Footprint: High - 30"w x 20.75"d
• Seat Dimensions: 27"w x 16.25"d	• Weight Capacity: 300 lb

Heavy Duty Transfer Bench
9075

- Converts easily for left or right hand transfers
- Adjustable height aluminum frame
- Latex-free rubber tips & suction tips



• Seat Height Adjustment: 16.75" - 21.75"	• Leg Footprint: High - 29.5"w x 20.75"d
• Seat Dimensions: 27"w - 16.25"d	• Weight Capacity: 500 lb

Padded Transfer Bench
9080

- Padded seat and back add comfort when bathing
- Converts easily for left or right hand transfers
- Adjustable height aluminum frame
- Latex-free rubber tips & suction tips



• Seat Height Adjustment: 17.75" - 21.75"	• Leg Footprint: High - 29.75" x 20.5"d
• Seat Dimensions: 27"w x 16"d	• Weight Capacity: 300 lb

Portable Transfer Bench
9072-R

- Great for travel and easy to store
- Tool-free assembly in minutes
- Converts easily for left or right hand transfers
- Latex-free rubber tips & suction cups



• Seat Height Adjustment: 19.25" - 23.25"	• Leg Footprint: High - 31.75"w x 22.5"d
• Seat Dimensions: 26.5"w x 16.25"d	• Weight Capacity: 300 lb

Padded Transfer Bench/Commode
9073

- Converts easily for left or right hand transfers
- Removable seat for easy cleaning
- Bucket and lid included
- Adjustable height aluminum frame
- Latex-free rubber tips & suction tips



• Seat Height Adjustment: 19" - 23.5"	• Leg Footprint: High - 36.5"w x 17.75"d
• Seat Dimensions: 26.75"w x 14.5"d	• Weight Capacity: 300 lb

Bedside Commode
8450

- Can be used beside the bed and over the toilet as a raised toilet seat with arms
- Adjustable height steel frame
- Bucket, lid, and splash guard included
- Toilet seat opening: 8.5" w x 9.5" d



• Seat Height: 17.5" - 23.5"	• Seat Dimensions: 13.25" w x 14.25" d
• Overall Dimensions: 24" w x 19.75" d	• Weight: 11.5 lb
• Width Between Arms: 17.75"	• Weight Capacity: 350 lb

Folding Commode
8700-R

- Can be used beside the bed and over the toilet as a raised toilet seat with arms
- Ideal for travel: folds and unfolds in seconds
- Adjustable height steel frame
- Bucket, lid, and splash guard included
- Toilet seat opening: 8" w x 10" d



• Seat Height: 16.25" - 22.25"	• Seat Dimensions: 13.75" w x 16.25" d
• Overall Dimensions: 22.5" w x 21.5" d	• Weight: 10 lb
• Width Between Arms: 18.5"	• Weight Capacity: 300 lb

Drop-Arm Commode
8900W

- Arms drop for easy patient transfer
- Bucket, lid, and splash guard included
- Available in Gray and Chrome
- Adjustable height steel frame
- Toilet seat opening: 7.75" w x 9.75" d



• Seat Height: 18.75" - 22.25"	• Seat Dimensions: 13.5" w x 16" d
• Overall Dimensions: 22.5" w x 20.75" - 22.75" d	• Weight: 19.6 lb
• Width Between Arms: 18"	• Weight Capacity: 300 lb

Padded Drop-Arm Commode
8901W

- Arms drop for easy patient transfer
- Padded seat and backrest
- Bucket, lid, and splash guard included
- Adjustable height steel frame
- Toilet seat opening: 8.25" w x 11" d



• Seat Height: 19.25" - 23.75"	• Seat Dimensions: 16.5" w x 16.5" d
• Overall Dimensions: 22.5" w x 21.25" - 23" d	• Weight: 20 lb
• Width Between Arms: 17.25"	• Weight Capacity: 300 lb

Heavy Duty Commode
8580

- Extra-large seat provides added support
- Bucket and lid included
- Adjustable height steel frame
- Toilet seat opening: 8" w x 10" d



• Seat Height: 19.75" - 23.5"	• Seat Dimensions: 23" w x 18.5" d
• Overall Dimensions: 32.25" w x 23.25" d	• Weight: 22.25 lb
• Width Between Arms: 28.5"	• Weight Capacity: 450 lb

Heavy Duty Commode
8582

- Extra-large seat provides added support
- Bucket, lid, and splash guard included
- Adjustable height steel frame
- Toilet Seat Opening: 8.25" w x 10.25" d



• Seat Height: 17.5" - 24.25"	• Seat Dimensions: 23.25" w x 18.75" d
• Overall Dimensions: 27" w x 23.75" d	• Weight: 16.6 lb
• Width Between Arms: 23.75"	• Weight Capacity: 500 lb

Heavy Duty Drop-Arm Commode

8583

- Arms drop for easy patient transfer
- Extra-large seat provides added support
- Bucket and lid included
- Adjustable height steel frame
- Toilet Seat Opening: 8.25"w x 10.25"d



- | | |
|-----------------------------------|--------------------------------------|
| • Seat Height: 19.75" - 23.75" | • Seat Dimensions: 23.25"w x 18.75"d |
| • Overall Dimensions: 30"w x 22"d | • Weight: 18.1 lb |
| • Width Between Arms: 25.75" | • Weight Capacity: 500 lb |

Drop-Arm Transport Chair Commode

8805

- Arms drop for easy patient transfer
- 4" swivel wheels with locking casters
- Padded seat fits over commode and converts into transport chair
- Bucket & lid included
- Removable, swing away footrests
- Bucket can be removed without lifting user or raising seat
- Toilet Seat Opening: 8.75"w x 10"d
- NOT FOR USE IN SHOWER



- | | |
|---------------------------------------------|-------------------------------------|
| • Seat Height: 18.75" | • Seat Dimensions: 17.75"w x 17.5"d |
| • Overall Dimensions: 21"w x 22.25"d x 38"h | • Weight: 22 lb |
| • Width Between Arms: 17.5" | • Weight Capacity: 250 lb |

Shower Chair/Commode

8800

- Waterfall armrests and padded seat
- 4" swivel wheels with locking casters
- Bucket and lid included
- Toilet seat opening: 8.5"w x 10"d



- | | |
|-----------------------------------------------|-------------------------------------|
| • Seat Height: 22.25" | • Seat Dimensions: 15.5"w x 16.75"d |
| • Overall Dimensions: 22.75"w x 22.5"d x 34"h | • Weight: 13.25 lb |
| • Width Between Arms: 19.5" | • Weight Capacity: 250 lb |

Deluxe Shower Chair/Commode

8801

- Waterfall armrests and padded seat
- 4" swivel wheels with locking casters
- Bucket and lid included
- Removable, swing away footrests
- Toilet seat opening: 8.5"w x 10"d



- | | |
|---------------------------------------------|--------------------------------|
| • Seat Height: 22.25" | • Seat Dimensions: 16"w x 16"d |
| • Overall Dimensions: 22"w x 24"d x 36.25"h | • Weight: 18.1 lb |
| • Width Between Arms: 18.75" | • Weight Capacity: 250 lb |

Commode Bucket Liner

PA-1082

- Medical grade bag with super absorbent pad that quickly transitions body fluids into a gel
- Eliminates any mess, odor or spills
- Provides convenience & safety for both user and caregiver
- Convenient drawstring closure for easier cleanup
- Pad absorbs 17 oz. of body fluids
- 20 bucket liners per box (comes retail packaged)



Commode Bucket Liner Display

PA-1082

- Cardboard display
- Includes 18 retail-packaged commode liner boxes (20 liners per box)
- Dimensions: 20"w x 14.25"d x 30"h



Scan here to watch our quick video about the Commode Bucket Liner Display and its easy assembly!



Toilet Safety Rails 8200-R



- One piece mounting flange for easy width adjustment
- Rails bolt securely to the toilet
- Waterfall armrest

- **Arm Height Adjustment:** 26.5" - 30.5"
- **Weight:** 3 lb
- **Width Between Arms:** 21" - 24"
- **Weight Capacity:** 250 lb

Toilet Support Rails 8201-R



- Provides added stability when sitting and rising
- Frame bolts directly to toilet
- Cushioned, easy to clean foam armrests
- Adjusts to fit standard and elongated toilets
- Easy to install

- **Arm Height:** 11.5"
- **Overall Dimensions:** 23.5"w x 16.25"d x 11.5"h
- **Width Between Arms:** 17.25" - 19.75"
- **Weight Capacity:** 300 lb

Toilet Safety Support Frame 8205-R



- Freestanding frame
- Slip-resistant base
- Padded foam handles
- For use with both standard and elongated toilets
- Includes magazine rack
- Easy, tool-free assembly

- **Overall Dimensions:** 29.75"w x 19"d x 25.5"h
- **Width Between Arms:** 20.5"
- **Handle Height:** 25.5"
- **Weight Capacity:** 300 lb

Toilet Seat Risers 8341-R & 8342-R



- Allows the original seat and lid to be used
- Bolts directly to the toilet for added safety
- Rubber pads help prevent seat from moving
- Adds 3.75" height to the toilet seat
- 8341-R (Elongated)
- 8342-R (Standard)

- **Overall Dimensions:** 14"w x 19.25"d x 3.75"h (8341-R)
- **Aperture Dimensions:** 10.25"w x 13"d (8341-R)
- 13.75"w x 17.25"d x 3.75"h (8342-R)
- 9.5"w x 11"d (8342-R)
- **Weight:** 2.4 lb (8341) & 2.3 lb (8342)
- **Weight Capacity:** 300 lb

Toilet Seat Risers 8343-R & 8344-R with Arms



- Allows the original seat and lid to be used
- Bolts directly to the toilet for added safety
- Rubber pads help prevent seat from moving
- Adds 3.75" height to the toilet seat
- Padded arms
- 8343-R (Elongated)
- 8344-R (Standard)

- **Overall Dimensions:** 21.25"w x 19.25"d x 9.75"h (8343-R)
- **Aperture Dimensions:** 10.25"w x 13"d (8343-R)
- 21"w x 17.25"d x 9.75"h (8344-R)
- 9.5"w x 11.25"d (8344-R)
- **Weight:** 3.7 lb (8343) & 3.55 lb (8344)
- **Weight Capacity:** 300 lb

Toilet Seat Risers (Hinged) 8345-R & 8346-R



- Hinged design makes cleaning easier and faster than before
- Allows the original seat and lid to be used
- Bolts directly to the toilet for added safety
- Adds 3.75" height to the toilet seat
- 8346-R (Elongated)
- 8345-R (Standard)

- **Overall Dimensions:** 13.75"w x 17"d x 3.75"h (8345-R)
- **Aperture Dimensions:** 9.5"w x 11.25"d (8345-R)
- 14.25"w x 19.25"d x 3.75"h (8346-R)
- 10.25"w x 14"d (8346-R)
- **Weight:** 2.7 lb (8345) & 3.05 lb (8346)
- **Weight Capacity:** 300 lb

5" Raised Toilet Seat

8340-R

- Adds 5" to the height of the toilet seat to aid in sitting and rising
- Non-locking
- Extra large seating surface



- | | |
|----------------------------------------------------|---------------------------|
| • Overall Dimensions:
17.5"w x 17.25"d x 6.25"h | • Weight: 3 lb |
| • Aperture Dimensions:
7.5"w x 9.5"d | • Weight Capacity: 300 lb |

5" Locking Raised Toilet Seat

8350-R

- Adds 5" to the height of the toilet seat to aid in sitting and rising
- Front locking mechanism for easy installation



- | | |
|----------------------------------------------------|---------------------------|
| • Overall Dimensions:
16.5"w x 17.25"d x 7.25"h | • Weight: 4.15 lb |
| • Aperture Dimensions:
7.75"w x 9"d | • Weight Capacity: 300 lb |

5" Locking Raised Toilet Seat

8351-R with Adjustable Width Detachable Arms

- Adds 5" to the height of the toilet seat to aid in sitting and rising
- Arms provide support for sitting and rising
- Front locking mechanism for easy installation
- Padded arms



- | | |
|---------------------------------------------------|---------------------------------------|
| • Overall Dimensions:
22.25"w x 17"d x 14.25"h | • Width Between Arms:
17.75" - 20" |
| • Aperture Dimensions:
7.25"w x 8.5"d | • Weight: 5.25 lb |
| | • Weight Capacity: 300 lb |

5" Locking Raised Toilet Seat

8352-R & 8353-R with Arms

- Adds 5" to the height of the toilet seat to aid in sitting and rising
- Front locking mechanism for easy installation
- Padded arms provide support for sitting and rising (8353-R only)



- | | |
|------------------------------------------------------|-------------------------------------------|
| • Overall Seat Dimensions:
20.5"w x 17.75"d x 7"h | • Width Between Arms: 18" |
| • Aperture Dimensions:
8"w x 9.5"d | • Weight: 4.6 lb (8352)
5.75 lb (8353) |
| | • Weight Capacity: 300 lb |

Padded Toilet Seat Risers

2629-R (2") & 2630-R (4")

- High density foam minimizes pressure on sensitive areas
- Attaches easily to toilet seat with four adjustable Velcro straps to reduce shifting
- Easy to clean vinyl surface
- Fits most toilet seats
- Dimensions:
2629-R: 14"w x 15.75"d x 2"h
2630-R: 14.5"w x 15.75"d x 4"h
- Weight Capacity: 250 lb



Easy Air Adjustable Toilet Seat Riser

2628A-R (2.25")

- Adjustable air foam combination
- Air release valve easily adjusts to your individual comfort
- Automatically re-inflates with no pump required
- Easily attaches to your toilet seat with loop straps
- Dimensions:
14"w x 15.5"d x 2.25"h
- Weight Capacity: 250 lb





PREVENT FALLS

This is one epidemic we can cure



IT'S NOT AN AGE THING

Fall related injuries in the bathroom are as common among those 18-25 as they are among those over 65. Let's make the bathroom safe for everyone!



Scan the QR Code to learn more about Fall Prevention and download NOVA's *Fall Prevention Guide*

Consult your physician before use of any medical products



Empowering Human Ability Since 1993

1-800-557-6682
novajoy.com

Authorized NOVA Dealer