

# MOBILITY


## TRANSPORT CHAIR WITH FLIP-UP DESK ARMS ITEM # 352BK (BLACK)



### THE FEATURES

- 1 Lightweight aluminum frame
- 2 Flip-up desk arms
- 3 Hand brakes & locking rear wheels
- 4 Non-skid foot plates & heel loops
- 5 Foldable
- 6 8" front wheels & 12" rear wheels
- 7 Anti-tippers

#### Available in:



We love the 352 Transport Chair because it's our most versatile transport chair, taking you the most places. It features a lightweight aluminum frame, flip-up desk arms and rugged 12" rear wheels. And the 352 just got even better, with a new, sleek color option in **BLACK**.

Item #	Weight Capacity	Overall Weight	Width btw Arms	Seat Height	Seat Dimensions	Overall Dimensions	Overall Folded Dimensions
352	300 lb	27 lb	19.75"	20"	18.5"w x 15.75"d	25.75"w x 31"d x 38"h	13"w x 29"d x 30.75"h

Los Angeles - Chicago

©2019 NOVA

CARE • QUALITY • INNOVATION



THE NOVA PROMISE

[novajoy.com](http://novajoy.com)

**CONTACT US!**  
We're here to help



[info@novajoy.com](mailto:info@novajoy.com)

1 800 557 6682