

Avenger

4 wheel

Seat Color Charcoal Vinyl

Body Colors Silver Re

> The Golden Avenger four-wheel, heavy duty scooter provides an impressive blend of rugged outdoor maneuverability while giving a smooth, comfortable ride over a variety of

nature's most challenging terrains. Featuring a 500 lb. weight capacity, the Avenger comes standard with NF22 batteries providing a range of 18 miles.





GA-541

Avenger 4 wheel **Specifications**

| * | |
|--|------------------------|
| Model | GA541 |
| Medicare Code | Not Classified |
| Weight Capacity | 500 lbs |
| Drive Wheels | Rear |
| Maximum Speed | 7 mph |
| Operating Range 1 | 18 miles |
| Ground Clearance to Center D | |
| Ground Clearance to Transaxle | |
| Turning Radius | 67" |
| Type Batteries | 2-NF22 |
| Colors | Silver, Red |
| Freewheel Mode | Yes |
| Electronic Speed Control | Yes |
| Electro-Mechanical Brakes | Yes |
| Charger | Off Board |
| Controller | PG 110 Amp |
| Dimensions: | |
| Length | 57.5" |
| Width | 24.5" |
| Height (ground to top | |
| of back of seat) | 42.5" - 45.5" |
| Height (ground to top | |
| of back of headrest) | 50" - 56" |
| Ground to Top of Seat | 25" - 28" |
| Ground to Top of Deck | 6.5" |
| | 18.75" - 21.75" |
| Number of Seat Height Posit | |
| Size of Increments | 1" |
| Front Axle to Rear Axle | 38.5" |
| Weight of Unit: (assembled | - |
| Front Half | 68 lbs |
| Rear Half | 68 lbs |
| Base Seet w/A mms | N/A |
| Seat w/Arms | 55 lbs |
| Batteries Tires: | 37 lbs ea Pneumatic |
| Front | 13" |
| Rear | 13" |
| Casters | N/A |
| | 3" |
| Rear Anti-Tippers Standard Seat: | |
| Slide Seat | Captain's Seat N/A |
| Back Height - no headrest | 21" - 24" |
| Back Height - no neadrest Back Height with headrest | 27.5" - 31.5" |
| Width x Depth | 20" X 17" |
| Color | Gray Vinyl |
| Warranty: | Gray Villyi |
| Frame | Lifetime |
| Drive Train | 2 Years |
| Electronics ² | 2 Years |
| LICCUOTICS | 2 16912 |

^{1.} Battery range will vary due to rider weight, drive surface and terrain. Foam-filled tires are optional. The operational mileage will be less if scooter is equipped with foam filled tires.
² Electronics warranty excludes batteries.



Easy to

Battery Meter

Easy to

Charger

Reach

Plug

Read

Features

- 500 lb weight capacity
- Speeds up to 7 mph
- Bright safety lights
- 22 NF batteries
- Up to 18 Mile Range
- 4-Way Flashers
- Great Warranty







Autolight

Key Switch

Rubber

Grips