



Shown with Power Elevating Seat


Standard Colors:



Crimson Red



Azure Blue



Galactic Grey

MEET OUR NEW COMPANION!

- NEW Tiller Design with Digital Dash & Console
- NEW Tiller Mounted USB Charging Port
- NEW High Back Stadium Seat with Quilted Diamond Stitch
- NEW Full Front & Rear Suspension
- NEW Larger Optional Battery Size for Increased Range
- NEW Galactic Grey, Azure Blue & Crimson Red Color Shrouds
- 3.5 Amp Super Charger for Maximum Battery Life

Companion 3-Wheel (Full-Size) GC340

^{GOLDEN} Companion

Golden's popular Companion 3-wheel scooter has a sporty new look and upgraded comfort features! This next-generation model has a well-designed, wider infinite adjustable tiller with a digital dashboard featuring a percentage battery life indicator and water-resistant buttons on the control panel, plus two convenient cup holders. It offers ample leg and foot room on the wide front deck and full front and rear suspension for a smooth ride. The luxurious high back stadium-style seat with a quilted diamond stitch pattern and extra padding elevates the level of comfort and style while on the go! The sliding swivel seat is 18" wide by 16" deep. Durable vibrant red, blue, and grey color shrouds add to the sportiness. Other features include easy disassembly, front and rear lighting package with side reflectors, and new wheels. Standard U1 batteries offer a 12 mile per-charge range. An optional 50AH battery offers extended range. The Companion is the ultimate value in a full-time, yet portable scooter.

Standard Features:

- Under Mat Troubleshooting Guide
- Delta Tiller with Auto Lock Mechanism
- Easy Reach Tiller Adjustment
- Adjustable Sliding Seat
- Adjustable LED Headlight
- U-1 Batteries



GOLDEN Companion

FULL-SIZE 3-WHEEL SCOOTER

GC340

Companion 3-Wheel Specifications

Model	GC340
Medicare Code	K0807
Weight Capacity	400 lbs
Drive Wheels	Rear
Maximum Speed	5 mph
Operating Range ¹	19 miles
Ground Clearance Under Deck	4.5"
Ground Clearance to Transaxle/Motor	2.5"
Turning Radius	45"
Type Batteries / Air Travel Safe	2-U1 / up to 50AH / Yes
Colors	Grey, Red, Blue
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Energy Efficient - Auto Turn Off (Sleep Mode)	Yes
Charger	Portable 3.5 Amp Super Charger Off Board
Controller	90 Amp Dynamic R Series
Dimensions:	
Length	46"
Width	24.5"
Ground to Top of Seat	23" - 25"
Ground to Top of Deck	6.5"
Deck to Top of Seat	16.5" - 18.5"
Number of Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	34.5"
Weight of Unit: (assembled)	207 lbs
Tires: Flat Free / Non Marking Solid	
Front	9"
Rear	10.75"
Front & Rear Anti-Tippers	2"
Standard Seat: Quilted High Back	
Slide Seat	Standard
Back Height (no headrest)	18"
Width x Depth	18" x 16"
Color	Black
Transportation Seat Lock	Yes, Self Locking
Seat Rotation	90 Degree either direction 360 Degrees
Optional Seat: Stadium Style	
Back Height (no headrest)	16"
Width x Depth	20" x 18"
Color	Black
Seat Options:	
Pneumatic Rise Seat	Yes
Power Elevating Seat - 6" with Lockout	Yes
Warranty:	
Frame	Lifetime
Drive Train	2 Years
Electronics ²	2 Years
Batteries	1 Year

¹ Battery range will vary due to rider weight, drive surface, terrain, and battery type.

² Electronics warranty excludes batteries.



Digital Dashboard & User-friendly auto-tiller adjustment



Infinite adjustable tiller with ergonomic design



Conveniently located drink holders



Front and rear suspension for a smooth ride



Disassembles easily!



Shown with Power Elevating Seat



Phone: 800-624-6374

Fax: 800-628-5165

www.goldentech.com

Find us on
Facebook

www.facebook.com/goldentechnologies

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

SHT-GC340-022624