

MOBILITY

TRANSPORT CHAIRS
WITH FLIP-UP DESK ARMS
ITEM # 332 / 349 / 352



THE IMPROVEMENTS

- 1 Seat belt length increased from 41" to 47"
- 2 Now comes with anti-tippers

THE FEATURES

- 1 Lightweight aluminum frame
- 2 Flip-up desk arms
- 3 Hand brakes (332 and 352 only) & locking rear wheels
- 4 Non-skid foot plates & heel loops
- 5 Foldable
- 6 8" front wheels & 12" rear wheels (332 and 352 only)

Available in:

 332B / 349B / 352B - Blue

 332R / 349R / 352R - Red



Anti-tippers

352

We love these Transport Chairs because they feature a lightweight aluminum frame, flip-up desk arms and locking rear wheels. And they just got even better, with an **extra 6" of seat belt length** and **anti-tippers** for a safer, more comfortable user experience.

NOVA Retail Price

\$399.95 (332) / \$299.95 (349) / \$359.95 (352)

Item #	Weight Capacity	Overall Weight	Width btw Arms	Seat Height	Seat Dimensions	Overall Dimensions	Overall Folded Dimensions
332	400 lb	27.75 lb	22"	20"	21.5"w x 15.5"d	38"h x 28.5"w x 30"d	38"h x 11"w x 30.5"d
349	300 lb	23 lb	19.25"	19"	18.5"w x 15.5"d	38.5"h x 25"w x 28.5"d	30.5"h x 12.5"w x 26.75"d
352	300 lb	27.5 lb	19.25"	19.75"	18.5"w x 15.5"d	38.25"h x 25.5"w x 31"d	30.5"h x 12.5"w x 29.5"d

Los Angeles - Chicago

©2018 NOVA

CARE • QUALITY • INNOVATION




THE NOVA PROMISE

novajoy.com

CONTACT US!
We're here to help

 info@novajoy.com

 1 800 557 6682