



ACCESS ALL LEVELS OF
YOUR HOME





ENJOY BRUNO'S SIGNATURE STAIRLIFT

Experience the ultra-smooth, quiet ride of a Bruno Elite curved stairlift. Bruno is the only manufacturer who handcrafts the rail to ensure a precision fit that artfully flows into every corner and maximizes open space on the steps.

With heavy-duty capacity, high-performance engineering and a Limited Lifetime Warranty, you gain the peace of mind that a Bruno Elite curved stairlift will allow you to stay safely in your home long into the future.



3 TOP CONSIDERATIONS

- ✓ Handcrafted in USA
- ✓ 400 lb. (181kg) Lift Capacity
- ✓ Limited Lifetime Warranty

MADE IN USA

Bruno engineers and manufactures our stairlifts in Wisconsin, USA. Each curved rail is carefully crafted by hand to ensure an optimal fit and ride. We take pride in our artisan touch!





CUSTOM FIT. HIGH-PERFORMANCE.

Handcrafted. Elegant, handmade rail delivers a precision fit and smooth ride while the covered gear rack gives clean, graceful appearance.

Easy-to-Use. Armrest control glides stairlift up and down. Two wireless call/send remotes are also included.

Adjustable. Change seat and footrest height and space between arms for optimal fit for your body.

Compact. Vertical, custom rail leaves maximum open space on steps. Fold the arms, seat and footrest when not in use to save even more space.

Safe. Seat swivels 90° at top of steps and locks in place for safe exit away from stairs. Seat belt and obstruction sensors further ensure safety.

Dependable. Battery-powered operation ensures performance, even in the event of a power outage.

[Learn More
bruno.com](https://bruno.com)

AVAILABLE OPTIONS

Add extended rail positions to park the stairlift safely away from the steps or out of view and an optional mid-park charging station for stairs with middle landings.

Get an extra-wide seat or larger footrest to comfortably accommodate larger stairlift riders.

Eliminate the need to lean over with a power folding footrest. Fold the seat and the footrest will automatically fold, too. Convert from manual to a power swivel seat and the chair will automatically rotate with the push of a button.

Select from four additional seat fabric color options to blend your stairlift perfectly into your home's decor.



Not all options are shown.



Verified Customer

"I love my Bruno stairlift. It works so smoothly... It is truly a blessing."

Sharon W.





TECH SPECS

Lift Capacity:
400 lb (181 kg)

Standard Seat Size:
18 in wide x 15 in deep (457 mm x 381 mm)

Larger Seat Option:
22.5 in wide x 16.5 in deep (572 mm x 419 mm)

Standard Footrest Size:
10 in long x 9 in wide (254 mm x 229 mm)

Larger Footrest Option:
12 in long x 11 in wide (305 mm x 279 mm)

Curved Rail:
Installs within 7.5 in (190 mm) of the wall

WARRANTY

Lifetime Warranty on Major Components -- for the lifetime of the Bruno Stairlift as long as the original purchaser owns the stairlift.

2-Year Warranty on all components after date of installation by an Authorized Bruno Dealer.

30-Days Labor Coverage after date of installation by an Authorized Bruno Dealer for any manufacturing defect.

ABOUT BRUNO

Veteran-Founded | Family-Owned | Made in USA

Bruno makes stairlifts and vertical platform lifts to help people better access homes, public buildings and businesses, and scooter lifts to assist people in easily transporting their mobility devices.

